

VISION ULTRA MID-WHEEL POWER WHEELCHAIR

P330



- 14" Knobby tires for optimum weight distribution provide superior traction.
- Patented articulation mechanism allows front wheels climbing up the obstacles. Plus the high traction of 14" tires, makes the vehicle drive over rough terrain easily.
- Dual in-line motors for enhanced efficiency, torque, range and performance.
- Full front and rear suspension design offers a comfortable ride.
- Mid-wheel drive design provides outstanding maneuverability in tight, compact spaces.
- Built to navigate everyday obstacles.
- Maximum stability with 4 seat post design.
- Swing back joystick mount is standard.

Driving Range	32km
Max Speed	7.2 km/h , 4.5 mph
Battery Type	12V/ 22NF x 2pcs
Motor	DC 24V , 200W(Rating) / 575W(Max)
Max Climbing Angle	10°
Turning Radius	53 cm / 21"
Ground Clearance	8 cm / 5"
Controller	PG nVR2-50A
Charger	5A off-board
Weight Capacity	205 kg / 450 lb
Total Weight	114 kg / 251 lb
Front wheel	8" x 2"PU tire
Drive Wheel	14"
Rear wheel	8" PU tire

Width 66cm	Seat Width 51cm	Seat Depth 46cm	Seat Height 56-58cm	Length 101cm	Footrest 43-48cm	Back Height 51cm	Height 127cm	Max Load 205kg	Product Weight 114kg

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za