

BASIK FOLDING SHOWER SEAT



FOLDING SHOWER SEAT



FOLDING SHOWER SEAT - WITH LEGS



DELUXE SHOWER SEAT



NAME	BASIK FOLDING SHOWER SEAT	FOLDING SHOWER MOULDED SEAT	FOLDING SHOWER SEAT – WITH LEGS	DELUXE FOLDING SHOWER SEAT
PRODUCT CODE	SHM5054	SMN6104	SMN5051	FSC
FRAME	POWDER COATED STEEL	POWDER COATED MILD STEEL	POWDER COATED MILD STEEL	STAINLESS STEEL
SEAT SIZE	37X23X2.7CM	490X430X480 (MM)	400X400 (MM)	460X390 (MM)
SEAT	BLOW MOULDED PLASTIC SEAT	ONE PIECE MOULDED PLASTIC SEAT	ONE PIECE MOULDED PLASTIC SEAT	ONCE PIECE SOLD PNENOLIC SEAT
WEIGHT		3.6KG, GROSS WEIGHT 4.2KG	4KG, GROSS WEIGHT 4.5KG	7KG, GROSS WEIGHT 7.5KG
LOAD CAPACITY	300LBS	80KG	100KG	90KG
PACKING	38X33X6CM	SEAT 37DX40D (CM) 50CM FROM WALL TO FRONT OF SEAT	490X600X150 (MM)	580X450X150 (MM)