

# BASIK COMMODE RANGE

## BASIK Taxi Commode

BASIK TAXI – 5" ALK694



- Chrome Steel Frame
- Drop Arm for Transfer
- 5" solid castors with brakes on rear castor
- Removable Footrest

|               |                    |                    |                     |                 |                     |                |                   |                    |                         |              |
|---------------|--------------------|--------------------|---------------------|-----------------|---------------------|----------------|-------------------|--------------------|-------------------------|--------------|
|               |                    |                    |                     |                 |                     |                |                   |                    |                         |              |
| Width<br>55cm | Seat Width<br>43cm | Seat Depth<br>43cm | Seat Height<br>51cm | Length<br>111cm | Back Height<br>47cm | Height<br>95cm | Max Load<br>110Kg | NW/GW KGS<br>14/16 | Packaging<br>62X57X55cm | PCS/CTN<br>1 |

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# BASIK COMMODE RANGE

## BASIK TAXI Aluminium Shower Commode

BASIK TAXI SHOWER – 5" CA613L



- Aluminium Powder Coated Frame
- Flip-up Armrest for transfer
- 5" Solid Castors with brakes on rear castor
- Removable Footrest

|               |                    |                    |                     |                 |                     |                |                   |                    |                         |              |
|---------------|--------------------|--------------------|---------------------|-----------------|---------------------|----------------|-------------------|--------------------|-------------------------|--------------|
|               |                    |                    |                     |                 |                     |                |                   |                    |                         |              |
| Width<br>55cm | Seat Width<br>46cm | Seat Depth<br>40cm | Seat Height<br>51cm | Length<br>111cm | Back Height<br>49cm | Height<br>93cm | Max Load<br>110Kg | NW/GW KGS<br>12/13 | Packaging<br>54X53X57cm | PCS/CTN<br>1 |

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# BASIK COMMODE RANGE

## BASIK Folding Commode

BASIK FOLDING COMMODE – ALK616L



- Aluminium Frame
- Folding Frame
- Height Adjustable

|                      |                           |                           |                               |                       |                            |                          |                          |                         |                                |                     |
|----------------------|---------------------------|---------------------------|-------------------------------|-----------------------|----------------------------|--------------------------|--------------------------|-------------------------|--------------------------------|---------------------|
|                      |                           |                           |                               |                       |                            |                          |                          |                         |                                |                     |
| <b>Width</b><br>56cm | <b>Seat Width</b><br>38cm | <b>Seat Depth</b><br>44cm | <b>Seat Height</b><br>45-55cm | <b>Length</b><br>52cm | <b>Back Height</b><br>35cm | <b>Height</b><br>78-88cm | <b>Max Load</b><br>100Kg | <b>NW/GW KGS</b><br>6/7 | <b>Packaging</b><br>50X22X73cm | <b>PCS/CTN</b><br>2 |

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# BASIK COMMODE RANGE

## BASIK Folding Commode – 3”

BASIK FOLDING COMMODE – 3” – ALK616L-3



- Aluminium Frame
- Folding Frame
- Height Adjustable
- 3” Solid Castors with brakes on rear castor

|               |                    |                    |                        |                |                     |                   |                   |                  |                         |              |
|---------------|--------------------|--------------------|------------------------|----------------|---------------------|-------------------|-------------------|------------------|-------------------------|--------------|
|               |                    |                    |                        |                |                     |                   |                   |                  |                         |              |
| Width<br>56cm | Seat Width<br>38cm | Seat Depth<br>44cm | Seat Height<br>44-55cm | Length<br>52cm | Back Height<br>35cm | Height<br>78-88cm | Max Load<br>100Kg | NW/GW KGS<br>6/7 | Packaging<br>50X22X73cm | PCS/CTN<br>2 |

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za

# BASIK COMMODE RANGE

## BASIK Standard Commode

BASIK STD COMMODE – COM001



- Steel Powder Coated Frame
- Fixed Frame
- Height Adjustable

|               |                    |                    |                        |                |                     |                   |                   |                    |                         |              |
|---------------|--------------------|--------------------|------------------------|----------------|---------------------|-------------------|-------------------|--------------------|-------------------------|--------------|
|               |                    |                    |                        |                |                     |                   |                   |                    |                         |              |
| Width<br>57cm | Seat Width<br>46cm | Seat Depth<br>44cm | Seat Height<br>45-55cm | Length<br>48cm | Back Height<br>35cm | Height<br>78-88cm | Max Load<br>115Kg | NW/GW KGS<br>12/15 | Packaging<br>50x22x73cm | PCS/CTN<br>4 |

# Medop<sup>cc</sup>

Medopark, 1 Manchester Road, Wadeville,  
P.O. Box 51124, Raedene, 2124  
enquiries@medop.co.za | www.medop.co.za  
JHB Tel: 011 827 5893/4/5  
CT Tel: 021 534 8628 | admin@medop.co.za