

BASIK COMMODE RANGE

BASIK Folding Commode

BASIK FOLDING COMMODE – ALK616L



- Aluminium Frame
- Folding Frame
- Height Adjustable

Width 56cm	Seat Width 38cm	Seat Depth 44cm	Seat Height 45-55cm	Length 52cm	Back Height 35cm	Height 78-88cm	Max Load 100Kg	NW/GW KGS 6/7	Packaging 50X22X73cm	PCS/CTN 2

Medop^{cc}

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za