

# HIP KIT



## Hip Kit Includes:

26" Reacher  
 Contoured Sponge  
 Sock & Stocking Aid  
 18" Flexible Plastic Shoehorn

### 86-0070 FabLife Hip Kit Includes:

- A. 50-1100 FabLife 26" Reacher
- B. 45-2381 FabLife Contour Sponge
- C. 86-0041 FabLife Deluxe Formed Sock Aid
- D. 86-0391 FabLife 18" Plastic Shoehorn



**A. 50-1100**  
**FabLife 26" Reacher**

- Extend your reach more than two feet
- Strong, molded rubber jaw with traction grooves for easy gripping
- Lightweight, durable frame
- Ergonomically designed handle
- Magnetic tip



**C. 86-0041**  
**FabLife Deluxe Formed Sock Aid**

- Designed to guide sock over heel
- Foam patch helps keep sock in place
- Rigid design
- 2 foam handles



**B. 45-2381**  
**FabLife Contour Sponge**

- 22½" L Plastic handle with foam sponges
- Has a hole in the end for hanging
- 5x3" hourglass shape sponge



**D. 86-0391**  
**FabLife 18" Plastic Shoehorn**

- Makes putting on shoes easier
- Flexible
- Convenient
- Economical