

QUADWALKER - GENESIS 3 EXTRA WIDE - HEAVY DUTY

Registered design

9201A

Oval tube frame construction



New Handle bars



Easy push button handle height adjustment



- Handgrips - ergonomic design
- Brakes - light touch lock brakes - suitable for arthritis sufferers

- Backrest - foldable for easy storage

- New Seat - extra wide double stitched seat - breathable nylon cover

- New Shopping bag has 75% more storage space

- Strengthened lightweight aluminium

- Wheels - large indoor/outdoor wheels fitted in a "U" shape fork

- Easy fold rounded backrest for comfort.
- For institutions, rental or domestic use.
- Suitable for indoor and outdoor use.
- Heavy duty construction - with a safe working limit of 200kg, Australian Standards Certified.
- New "light touch" lock brake system.
- The unit can be folded without removing the storage bag.
- New nylon bag prevents visibility of valuable items. e.g. Wallet, Hand Bag, Keys

Rear frame widened for maximum walking space



Shopping bag has 75% more storage space



New bigger and wider wheels suitable for outdoor/indoor use



9201A - Heavy duty - extra wide

overall height	overall length	overall width	seat depth	seating width	seat height	folded width	castor size	unit weight	weight capacity
min 87cm max 99cm	70.5cm	71cm	35cm	56cm	55cm	24.5cm	20cm	9.2kg	200kg

Medop cc

Medopark, 1 Manchester Road, Wadeville,
P.O. Box 51124, Raedene, 2124
enquiries@medop.co.za | www.medop.co.za
JHB Tel: 011 827 5893/4/5
CT Tel: 021 534 8628 | admin@medop.co.za